


Use a jug from soap or sanitizer (heavier than average), <u>not</u> pesticide. Replace hose or jug as needed. Can use without clamp device if jug kept below tube tip until inserted in calf's esophagus. Can avoid cutout in jug if wide-mouth jug is available.





FOR MORE INFORMATION, CALL OR WRITE MILK PRODUCTS, LLC P.O. Box 150, Chilton, Wisconsin 53014, (920) 849-2348, Fax (920) 849-9014 www.milkproductsinc.com

FL/T.33e © 2009 Milk Products, LLC